Annolid: an instance segmentation-based multiple-animal tracking and behavior analysis package


Chen Yang, Matthew Einhorn, Thomas A. Cleland Dept. of Psychology, Cornell University, Ithaca, NY 14853




Introduction

Background: Animal tracking with computer vision

- 1. DeepLabCut and SLEAP: animal pose estimation
- 2. SimBA (behavior recognition from pose)
- 3. Annolid: instance segmentation







https://github.com/DeepLabCut/DeepLabCut https://github.com/murthylab/sleap

Why do we need instance segmentation?

Annolid treats keypoints-based pose estimation as a special case of instance


segmentation.

Instance mask area and perimeter can be tracked, as well as keypoints.

Object Investigation




Area of ROIs



Annolid uses deep learning for animal tracking and behavior analysis

<u>Annolid</u> = <u>Annotation</u> + <u>Annelid</u> (segmentation)




Original video

Mask
(Mask R-CNN or YOLACT)

Behavior identification


Annolid workflow


- Labeling of frames (annotation)
- COCO formatting
- Training and inference (local or Colab)
- Post-processing and analysis



How do we label animals and behaviors?

Polygons & keypoints(e.g. Vole_1, nose, tail_base, rearing.....) ~100 labeled frames






Tracking and re-identification

To link instances across frames, we treat each instance as its own class across frames.

To track multiple animals, label each animal as a separate instance.

To generalize across animals or videos, label multiple animals as examples of the same instance.



Animal tracking

Keypoint tracking of body parts

Automated behavior identification (bounding box)



Multiple animal tracking, including periods of partial body occlusion

Whole-body masking

Automated identification of interactions (bounding box)



Animal tracking

Masking and automatic scoring of *lone* animals and *huddles* of multiple animals


Confidence of identification reported



Animal and object tracking, including periods of occlusion


Tracked objects automatically associated with user-defined zones

Robustness to noisy background



Identification of freezing behavior (e.g., from fear conditioning)


Reporting of *motion*score based on optical
flow measurements
applied selectively to the
body mask



Multiple animal tracking on cryptic background



Multiple animal tracking with a large field of view



Accessibility and efficiency

- Options for training on Google Colab (as well as on a local workstation)
- Fast training with quality- and speed-optimized options


Model training

- 200 labeled images
- < 2 hours for 3000 iterations on Colab
- 30 min on NVidia 1080Ti

Inference (applying trained model to behavior videos)

- Mask R-CNN: ~7 FPS
- YOLACT: ~30 FPS

Transfer learning based on existing models trained on COCO dataset improves performance



Capacities for autolabeling and human-in-the-loop iterative model training

Conclusion and future directions

- Annolid is based on *instance segmentation* models, which facilitate the tracking of multiple animals along with flexible state identification (e.g., behavior classification, urine deposition, interactions among objects)
- Self-supervised, weakly-supervised, and unsupervised training options
- Incorporation of optical flow mechanics for improved performance
- Improved labeling efficiency via autolabeling and iterative model training

Annolid and Glitter2 are both freely available:

https://cplab.science/annolid https://cplab.science/annolid-tutorial-det2 https://cplab.science/glitter2



Thanks for listening!